



# *One Step Dementia Risk Test Kit*

(Colloidal Gold Lateral Immunochromatography)



**Hunan Qankorey Biotechnology Co., Ltd.**



## Features



Accurate



Non-Invasive



Easy



Fast 10 Minutes



Clinically Approved



Globally First



Urine-based amyloid- $\beta$ (A $\beta$ ) biomarker testing for risk of Alzheimer’s Disease

## Specifications

Product Name	One Step Dementia Risk Test Kit (Colloidal Gold Lateral Immunochromatography)
Certificate	CE, NMPA, MDA
Packing	1pcs/box, 25pcs/box, 50pcs/box
Usage	This product is an <i>in vitro</i> auxiliary screening device for risk alert of dementia (including but not limited to Alzheimer's disease). It may be used to evaluate brain health status and assist doctors to decide whether further tests are required.



# Clinical Performance

75.7%

Sensitivity

91.03%

Specificity

0.67

Kappa Value

Diagnostic group	Qty	Aβ-PET Scan positive	Urine Aβ detection positive	Aβ-PET Scan positive coincidence rate	Urine Aβ detection positive coincidence rate
AD group	80	75	64	93.75%	80.00%
Other dementia	6	1	4	16.7%	66.7%
MCI	21	10	13	47.62%	61.90%
CN	301	/	27	/	8.97%

Reference Principle				In total
	Case Group( + )	Comparison Group( - )		
Urine Aβ Detection	Positive (+)	81	27	108
Urine Aβ Detection	Negative ( - )	26	274	300
In total		107	301	408

Sensitivity=A/(A+C)=81/107x100%=75.7%

Specificity =D/(B+D) =274/301x100%=91.03%

Actual Coincidence Rate  $P_o = \frac{(A+D)}{(A+B+C+D)} \times 100\% = \frac{(81+274)}{408} \times 100\% = 87\%$

Theoretical Coincidence Rate  $P_e = \frac{(A+B)*(A+C)+(C+D)*(B+D)}{(A+B+C+D)^2} \times 100\% = \frac{108 * 107 + 300 * 301}{408^2} \times 100\% \approx 61.2\%$

Kappa Value  $K = \frac{P_o - P_e}{1 - P_e} = \frac{0.87 - 0.612}{1 - 0.612} = 0.67$

This clinical verification was completed independently by Sichuan Provincial People's Hospital in 2023. Total 408 participants were enrolled and classified as the cognitively normal (CN) group, the mild cognitive impairment (MCI) group, the AD group and the other dementia group. Meanwhile, Aβ-PET Scan was performed for comparison in clinically defined patients. The results support our One Step Dementia Risk Test Kit as an early screening device for MCI and dementia, and warrant the further longitudinal studies.



## Instructions

### 1 Drinking

Please drink 200-300mL water 30-60 min before test.



### 3 Testing

Draw some urine with the supplied straw, and slowly dispense 4 drops, one by one, to the sample hole of the test strip.



### 2 Sampling

Collect midstream urine sample with the disposable cup supplied.



### 4 Result

Place the strip on a table, wait for 10 min, then take a picture.




### 5 Positive


If positive, repeat 1 or 2 times in weeks. If the positivity continues, it is recommended to see a doctor for further check.



### 6 Negative & Invalid



If negative, monitor again in every 6-12 months.



Please contact customer service.

## Interpretation of test results

**Positive:** The concentration of amyloid- $\beta$  protein(A $\beta$ ) is above the normal range, indicating the high risk of AD progression. It is recommended to check the potential risk factors accordingly and adjust life style, e.g.(Reducing alcohol and smoke intake, getting quality sleep etc.)

**Negative:** The concentration of amyloid- $\beta$  protein(A $\beta$ ) is in the normal range. It is recommended to maintain the existing lifestyle and test regularly every year.



**Hunan Qankorey Biotechnology Co., Ltd.**

**Website:** <http://en.qankorey.com>

**Contact:** [lisa@qankorey.com](mailto:lisa@qankorey.com)

**WhatsApp:** +86 188 7423 2226